Curapuls 970

Bedieningshandleiding
Operating Instructions
Gebrauchsanweisung
Mode d'emploi
Instrucciones de manejo

Copyright:



Enraf-Nonius B.V.
P.O. Box 810
NL-2600 AV DELFT
Netherlands
Tel.: +31 (0)15 - 26 98 400
Fax: + 31(0)15 - 25 61 686
www.enraf-nonius.com

Article number 1419.751_42
Feb 2002

Curapuls 970

Bedieningshandleiding
Operating Instructions
Gebrauchsanweisung
Mode d'emploi
Instrucciones de manejo

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INTRODUCTION

The Curapuls 970 is a shortwave unit that enables HF energy to be applied selectively deep in the tissues and joints without unnecessarily affecting the surrounding area. This is effected through capacitive or inductive electrodes. The capacitive electrodes have a larger surface effect, the induction electrodes a larger effect in the deeper tissues. The equipment is only meant to be used by competent personnel in physiotherapy, rehabilitation or adjacent areas.

The Curapuls 970 offers two therapy possibilities: continuous and pulsed.

Continuous shortwave therapy

In continuous shortwave therapy the subjective feeling of warmth plays the major role. The output is continuous and can be set to a maximum of 400 W. The automatic tuning ensures that the selected intensity is maintained throughout the treatment. This technique increases the circulation and stimulates metabolic processes.

Pulsed shortwave therapy

Pulsed shortwave therapy is an athermic technique with enhanced biological effects. The maximum output is 1000 W (peak power), the mean output is much lower: between 6 and 80 W, depending on the selected pulse repetition frequency (which can be set in 10 steps from 15 - 200 Hz).

Furthermore, acute conditions where heat treatment is contraindicated can very well be treated with shortwave therapy. Even the presence of metal implants is not an absolute contraindication.

Some specific indications:

- Post-traumatic conditions.
- Post-operative conditions.
- Acute inflammations.
- Acute sport injuries such as hematomas and sprains.

See further indications and contra-indications on page 21.

Other characteristics of the Curapuls 970:

- Rectangular pulses.
- Fast, solid-state automatic tuning.
- Provision for HF magnetic field therapy.
- Multi-directional articulated electrode arms with single-knob fixation.
- Large castors with safety locks.
- Wide range of accessories (see our catalogue).

USER'S RESPONSIBILITIES**General**

The Curapuls 970 is a modern, safe and easy-to-handle apparatus, which has been developed, produced, calibrated and packed with the greatest care. The Curapuls 970 meets the requirements of the international standard for the safety of medical electrical equipment IEC 60601-1 (General standard), IEC 60601-2-3 (Standard for shortwave therapy equipment) and IEC 60601-1-2 (Standard for EMC). According to the IEC 60601-1-2 (1993) directive, clause 36.201.1.8, class A equipment and/or systems (CISPR 11 Classification) are allowed in domestic establishments when used under the jurisdiction of a health care professional.

An Enraf-Nonius test certificate is supplied with the unit. This certificate gives important measurements with respect to safety and specifications.

Enraf-Nonius cannot be held responsible for any discomfort to the operator or the patient due to faulty diagnosis, misuse or mishandling of the equipment and/or accessories, misunderstanding of the operating instructions or maintenance errors.

Service and maintenance

To meet the legal regulations of the countries that have implemented the IEC 60601-1 standard, it is necessary to have the Curapuls 970 checked and safety tested once each year. We advise keeping your equipment under service/maintenance at your local authorized service office. It is also recommended that a record be kept of all activities related to service and maintenance. In some countries this is even obligatory. The check up and/or technical maintenance must be carried out conform the procedure described in the service manual of the unit.

Unqualified personnel should not open the equipment for any reason whatsoever.

Enraf-Nonius will not be held responsible for damage to equipment or injuries to persons operating or being treated with the equipment that result from service or maintenance performed by companies or persons not authorized by Enraf-Nonius.

Combination therapy

Since the Curapuls 970 radiates an electromagnetic field that can affect electronic circuits inside equipment in the direct vicinity, you are strongly advised against any sort of combination in treatment whereby a second (electromedical) apparatus is used.

Installation

Good air circulation is essential to prevent internal heat build-up.

Do not install the unit in a location near heat sources such as radiators or air ducts, or in a place subject to direct sunlight, dust, excessive mechanical vibration or shock. It is not recommended to use the Curapuls 970 immediately after extreme temperature fluctuations.

The Curapuls 970 is not intended for use in so-called "wet rooms" (hydrotherapy rooms).

When moving or transporting the unit, never pull at the electrodes or the electrode arm, to avoid the risk of tilting the Curapuls 970.

For safety reasons, the brakes on the wheels of the unit should be locked before the beginning of a treatment.

Should any liquid or foreign objects fall into the cabinet, turn off the power, unplug the unit and have it checked by your nearest authorized service dealer before any further operation.

Connection

The Curapuls 970 is a class I, type BF apparatus according to IEC 60601-1. This means that the Curapuls 970 must be connected to a grounded power supply.

Prior to connection to the mains, check that the voltage and frequency stated on the Curapuls 970 correspond with the available power supply rating. Check that the available power supply complies with local requirements regarding medical rooms (electrical codes).

PRODUCT LIABILITY

A law on Product Liability has become effective in many countries. This Product Liability law implies, amongst other things, that once a period of 10 years has elapsed after a product has been brought into circulation, the manufacturer can no longer be held responsible for possible shortages of the product.

To the maximum extent permitted by applicable law, in no event will Enraf-Nonius or its suppliers or resellers be liable for any indirect, special, incidental or consequential damages arising from the use of or inability to use the product, including, without limitation, damages for loss of goodwill, work and productivity, computer failure or malfunction, or any and all other commercial damages or losses, even if advised of the possibility thereof, and regardless of the legal or equitable theory (contract, tort or otherwise) upon which the claim is based. In any case,

Enraf-Nonius's entire liability under any provision of this agreement shall not exceed in the aggregate the sum of the fees paid for this product and fees for support of the product received by Enraf-Nonius under a separate support agreement (if any), with the exception of death or personal injury caused by the negligence of Enraf-Nonius to the extent applicable law prohibits the limitation of damages in such cases. Enraf-Nonius can not be held liable for any consequence resulting from incorrect information provided by its personnel, or errors incorporated in this manual and / or other accompanying documentation (including commercial documentation). The opposing party (product's user or its representative) shall disclaim Enraf-Nonius from all claims arising from third parties, whatever nature or whatever relationship to the opposing party.

REMARKS CONCERNING TREATMENT

General

We strongly advise that you regularly consult the therapy manual (optionally available) with the Curapuls 970. The chapters regarding indications and contraindications are especially important with respect to safe and effective treatment.

Also pay attention to the following:

- From considerations of safety, only original electrode cables and electrodes from Enraf-Nonius should be used.
- Electronic hearing-aids must be taken off unless damage or malfunctions are not to be expected.
- Patients under treatment must not be able to touch surrounding metal objects such as metal furniture, radiators or window-frames.
- The electrode cables must run free from surrounding objects and must maintain sufficient distance with respect to each other as well as with respect to the treated patient.
- Use cable-clamps to lead the cable along the electrode arm.
- The patient must be in sight of the therapist at all times.

Furniture and clothing

Metal parts in furniture might cause intensity-concentrations in the electromagnetic field, therefore only use metal-free wooden furniture.

For all treatments with shortwave it is advisable to undress the patient, and to apply a towel between the skin and the electrode (against perspiration). There are several advantages in treating naked patients: in the first place, it is the easiest way to ensure that there are no metals in the electro-magnetic field (metals in the clothing or in ornaments under the clothing), and secondly it prevents undesired concentrations of energy due to perspiration or the wearing of for example nylon or leather clothing.

Implanted electronic devices

Any patient with an implanted electronic device (for example a cardiac pacemaker) should not be subjected to any treatment with shortwave equipment and should not remain in the vicinity of a machine in operation.

Implanted metals

Metal concentrates electromagnetic energy. To prevent any undesirable concentrations around a metal implant and the resulting dangers (burns), continuous shortwave therapy should basically not be used. Consult also the therapy manual 'Pulsed and Continuous Shortwave Therapy' (chapter on Relative Contraindications).

Electromagnetic interference

The Curapuls 970 produces an electromagnetic field that is designed to generate thermal and biological effects in tissues and joints of the human body. In free air, the strength of this electromagnetic field generated from the electrodes, cables and treated patient is reduced quadratically with the distance. This electromagnetic field can affect electronic circuits of equipment positioned in the vicinity.

This electromagnetic interference is in almost all cases the cause of so-called direct radiation on the statutory permitted frequency of 27.12 MHz. Despite the fact that the affected equipment is not tuned to this frequency, detection appears in several circuits due to the strong radiation, causing the modulation of the electromagnetic wave to become audible and/or visible (through the neon checklight). With low- and medium-frequency physiotherapy equipment this modulation might become perceptible to the patient as well.

The susceptibility to electromagnetic interference of equipment from the various manufacturers can be very different. This is due to the fact that some manufacturers are still shielding their equipment insufficiently against external electromagnetic fields. Should interference arise, we advise consulting the manufacturer concerned.

We recommend paying attention to the following items in order to minimize the chance of electromagnetic interference.

- The larger the distance between the short wave unit and other equipment, the smaller the electromagnetic interference. Place your shortwave unit and the patient under treatment, at least 2 metres from surrounding equipment and/or oblong metal objects.
- The capacitive electrodes must be placed at a distance of less than 8 mm from the skin. In this way the leakage current is minimal.
- Make sure that no oblong conductive (metal)

objects are present in the direct vicinity because they can operate as secondary antennae (especially objects with a length of approximately 5.5 metres).

- Use a separate mains group (phase) for the shortwave apparatus. Keep mains leads wide apart.
- Never coil up a mains lead for shortening when a shortwave apparatus is operating in the vicinity, or use extension cords.
- Take care that the mains lead of the shortwave apparatus does not come in the vicinity of the patient and/or electrode leads. Either spread the lead out on the floor or put excess lead underneath the shortwave unit.
- Use only a completely metal-free treatment chair with footrest or a wooden couch when treating with shortwave diathermy equipment.
- Use an optimal treatment method. Make sure that your patient is comfortable and relaxed and that electrodes and cables are properly fixed.
- HF interference will be considerably reduced when the Circuplode magnetic field electrode is used. When the Curapuls 970 is placed in a **Faraday cage** HF interference will be completely eliminated.
- The intensity should only be adjusted when a patient is to be treated and when the electrodes have been placed.

(See fold out page at front)

[1] Time display

The treatment time is shown in half minutes. The flashing decimal point indicates that short wave energy is applied and that the timer is running.

[2] Intensity display

The bottom field lights when a minimum intensity is set. The higher the set intensity, the more fields will light.

[3] Indicator lamp for tuning

As soon as the apparatus is tuned, this lamp switches off. During treatment this lamp may briefly light if the patient moves.

[4] Intensity control

With this control the required intensity can be set. The maximum intensity ranges which may be set for the various electrodes are indicated above the control.

CPL When the Circuplode is used for applying **continuous** shortwave therapy, the intensity set should not exceed position 6.

FPL When the Flexiplode is used, the intensity set should not exceed position 7. This applies to continuous as well as to pulsed shortwave therapy.

CPL When the Circuplode is used for applying **pulsed** shortwave therapy, the intensity set should not exceed position 8.

First set intensity. Then turn the intensity control completely to the left until a 'click' is heard ('0'-position); from here the required intensity can be set by turning to the right.

[5] and [6] Time setting switch

The treatment time can be set from 0 to 30 minutes. Set time: press the "+"-key.

Reduce time: press the "-"-key.

By pressing the "+"-key one time, the treatment time can be set from 0 to 30 min. in 1 step. The time starts to run back to '0' as soon as the intensity is set. After expiration of the treatment time the HF generator is automatically switched off and a warning buzzer sounds. Treatment can only be started after setting treatment time [4].

[7] Selector switch for therapy form

Position "I" for continuous shortwave therapy.

Other positions for pulsed shortwave therapy.

The following pulse repetition frequencies can be set: 15 - 20 - 26 - 35 - 46 - 62 - 82 - 110 - 150 - 200 Hz.

[8] Mains switch

Switch on: press the switch to (I).

Switch off: press the switch to (0).

When switched on, treatment time display [1] and indicator lamp for tuning [3] should light.

[9] Attachment points for electrode arms

By means of the supplied tool the electrode arms are mounted in these holes (see 'Electrode arms').

[10] Earth socket for adaptor

Via this socket the Circuplode can be earthed.

[11] Sockets for electrode cables

All the electrodes are connected to these sockets, with or without adaptor (see 'Connection of electrodes').

Connection of cables or accessories other than those prescribed by Enraf-Nonius can adversely affect the safety of the patient and the functioning of the unit, and is therefore not permitted.

At the back of the unit you will find:

- Automatic resettable cut-out
This switch is normally in position 'I'. If an overload occurs, the cut-out switches off the HF-output. After approx. 2 minutes the apparatus can be switched on again by resetting the switch to the position 'I'. Set the clock to '0' before switching on again.
- Ventilation grate
This grate is necessary for cooling the apparatus, and must **not** be covered.
- Storage compartment
This compartment is intended to hold for instance the operating instructions and the test certificate.
- Warning plate
This warning plate indicates that dangerous high voltages are used in the apparatus.
- Type plate
This plate indicates data such as type, series and serial numbers (for warranty, service, etc.). Power supply data such as mains voltage and current consumption are also indicated.
- Connection for potential-equalization cable
If prescribed by the local authorities, a potential-equalization cable can be supplied.

- Approvals stickers
In this place, safety stickers can be affixed.
- Hook for mains cable
For hanging up the mains cable when the apparatus is not in use.

If required by the local high-frequency interference regulations, the apparatus can be provided with a **shielded mains cable** and an extra **interference filter**.

The metal fixation plate at the broad end of the electrode arm(s) is to be connected to the back of the Curapuls by means of the supplied hollow-head screw wrench. You can connect the arm with the knobs either towards the outside, or towards the unit.

The fixation knobs can be adjusted by hand. We advise you to adjust so tight that the arm can still be adjusted, without having to untighten the knobs.

In the top and bottom of the broad side of the electrode arm are holes in which you can click the (supplied) cable guides. If desired, you can also apply the blind plugs.

Should the purple-ish shell be removed from the small end of the arm, replace it and fix it by means of the fixation knob at the small end of the arm, before mounting an applicator. You cannot mount an applicator, when the purple-ish shell is not in its place.

CONNECTION OF ELECTRODES

For safety reasons only original electrodes and electrode cables must be used.

Capacitive electrodes*Disc electrodes (Schliephake).*

This essential metal electrode can be adjusted continuously by pulling out the dark grey rear part of the electrode: indicating grooves become visible (the adjustment is 5 mm per groove).

With these electrodes a towel should be placed between the skin and the electrode, after clearing the area to be treated.

Rubber plate electrodes.

Connect these via the sockets for electrode cables. With these electrodes the 2 cm thick felt spacers (see catalogue) should be used between the skin and the electrode.

Induction electrodes

With these electrodes a towel should be placed between the skin and the electrode, after clearing the area to be treated.

Flexiplode

Connect this electrode to the sockets for electrode cables.

Circuplode

Connect this electrode to the sockets for electrode cables via the Circuplode-adaptor.

Two way plug

With this plug, two capacitive electrodes can be combined with a rubber electrode.

OPERATING SEQUENCE

- a. Insert the mains plug in the earthed wall socket.
- b. Check that the resettable cut-out is in the upper position ('I').
- c. Connect the required electrode(s) (see our therapy manual).
- d. Switch on the unit [8]. The time display [1] and indicator lamp for tuning [3] must light.
- e. Place the electrode arms and the electrodes in treatment position.
- f. Check the connection of electrodes and cables (**ensure that there is sufficient distance between the various cables!**).
- g. Set the selector switch for therapy form [7] to the required position (continuous or pulsed). If the pulsed therapy is chosen, set the required pulse repetition frequency.
- h. Set treatment time [5,6] (this will start the ventilator).
- i. Turn the intensity control [4] completely to the left until a 'click' is heard ('0' position). The point of the time display [1] starts flickering and the bottom field of the intensity display [2] lights. Turn the intensity control [4] slowly to the right: after passing position '1' the indicator lamp for tuning [6] will go out. The unit is now tuned. Set the required intensity by turning the control [8] further to the right.
- j. When treatment time is ended, the intensity is automatically switched off and a warning signal sounds. The display [2] extinguishes.

Remark: always set intensity out of the '0' position before starting the treatment.

Consult our therapy manual for applications.

Operational errors

If you switch over from continuous to pulsed shortwave **during treatment:**

The automatic security cut-out switches off the HF generator. Lamp [3] lights up.

If you switch over from pulsed to continuous shortwave **during treatment:**

The pulse intensity is reduced to a safe minimal value.

In both cases reset as follows:

- First set treatment time back to '0' [5,6].
- Select required therapy form [7] (continuous or pulsed).

Repeat operating sequence as from item h.

Indications

Specific indications for pulsed shortwave therapy are:

Pain with:

- Arthrosis
- Neuralgias
- Neuritis
- Vasomotoric head ache
- Post traumatic conditions

Hypertonia with:

- Arthrosis
- Neuralgias
- Reflexhypertonia

Post-traumatic conditions, for example

- A sprain or dislocation
- Contusion
- Rupture
- Fracture
- Hematoma
- A graze

Post-operative conditions, for example

- jaw-, foot-, and hip operations

Inflammation, for example

- chronic osteitis
- bursitis (possibly with calcification)
- sinusitis

Peripheral blood circulation disorders

- Decutius
- Oedema

Absolute contra-indications

- Rheumatoid arthritis
- Malignant tumors
- Pacemakers
- Pregnancy
- Tuberculosis
- Fever / acute infections

Relative contra-indications

- Implanted metals
- Serious arterial and venous circulation disorders such as atherosclerosis, thrombosis etc.
- Heart disorders
- Acute infectious diseases, acute inflammation
- Thermic sensibility disorders
- Neuropathies

Unproved, largely traditional contra-indications

- Osteoporosis
- Rapidly deviding tissue
- Hemophilia
- Use of anticoagulant drugs

FAULT CONDITIONS

After switching on the Curapuls 970, the intensity display and indicator lamp for tuning does not light.

Check that there is power to the wall socket.

Indicator lamp for tuning does not go out, but intensity display lights.

Check electrodes and cables, and their connections.

Intensity display does not light, but indicator lamp for tuning lights.

Check that the automatic cut-out is in position '□'.

MAINTENANCE BY THE USER**General**

Before cleaning and other maintenance of the Curapuls 970, the unit must first be switched off, and disconnected from the mains supply.

Opening of the equipment and attempts at repair by unauthorized persons can adversely affect the safety and are therefore not permitted.

Apparatus

The apparatus may be cleaned with a damp cloth using a non-abrasive, liquid household cleanser. Also the electrode cables, electrode arms, plate electrodes and induction electrodes may be cleaned in a similar way.

Flexible rubber electrodes

Clean the flexible rubber electrodes, preferably daily, with normal lukewarm water. If necessary, a household soap may be used. To avoid corrosion of the cable junction, do not completely submerge the electrodes. Avoid contact of electrodes with grease or oil.

Before using the flexible rubber electrodes after cleaning, make sure that both surfaces are fully dry.

The electrodes may be disinfected with a cloth dampened with 70% alcohol. Be aware that the black of the rubber pad electrodes may stain.

Electrode cables

Check all cables and electrodes regularly for damage (e.g. minute cracks) and bad contacts. Please realize that any poor contact might cause a local overheating effect. Pay special attention to possible cracks in the insulation of the cables. Always replace damaged cables or electrodes immediately. We advise keeping an extra set of frequently used accessories in stock.

End life span apparatus and accessories

Your Curapuls 970 and accessories contains materials which can be recycled and/or are noxious for the environment. At the end of the duration of life, specialised concerns can take apart these articles and sort out the noxious materials and materials for recycling. By doing so you contribute to a better environment.

Please ensure that you are well informed of the local rules and regulations regarding to the removal of apparatus and accessories.

Generator frequency	: 27.12 MHz, $\pm 0.6\%$
Output power	: continuous HF max. 400 W */** pulsed HF max. 1000 W (peak) */**
Pulse duration	: ca. 400 μ s
Pulse repetition frequency	: 15 - 200 Hz adjustable in 10 steps
Timer	: 0 - 30 min.; with acoustic signal and automatic switch-off

** Models for Canada and the USA have different specifications. Please consult your supplier.

* Measured with 70 Ohm standard phantom (art.nr. 1462.471).

TECHNICAL DATA

Mains supply	: 220 - 240 V / 50 Hz (art.nr. 1419.942) : 240 V / 50 Hz (art.nr. 1419.962) : 115 V / 50 Hz (art.nr. 1419.982)
Mains voltage variation	: max. $\pm 10\%$
Current consumption	: max. 6 A (at 230 V)
Medical device classification	: IIa (according to Medical Device Directive (93/42/EEC))
Leakage current	: typically 350 μ A (IEC requirement $\leq 500 \mu$ A)
Safety class	: I, type BF according to IEC 60601-1
Dimensions	: 56 x 43 x 92.5 cm (w x d x h) without arms
Length electrode arms	: 93 cm
Weight	: ca. 80 kg

Environment conditions for transport and storage

Environment temperature	: -10° till +50° C
Relative humidity	: 10 till 95 %
Atmospheric pressure	: 500 till 1060 hPa

Environments conditions for normal use

Environment temperature	: 10° till 40° C
Relative humidity	: 10 till 90 %
Atmospheric pressure	: 500 till 1060 hPa

Classification

Medical class IIa

This equipment complies with all requirements of the Medical Device Directive (93/42/EEC).

International safety standard IEC 60601-1

Non-ionising radiation

I

Safety class I

The unit must be connected to a grounded wall socket.

Type BF (including Circuplodes)

The leakage current and the ground connection of the equipment meets the requirements of IEC 60601-1 and IEC 60601-2-3.

English

For the ordering data of the Curapuls 970, standard accessories and additional accessories we refer to the Catalogue for Physiotherapy.

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