Sonopuls 490 / 491

Bedieningshandleiding
Operating Instructions
Gebrauchsanweisung
Mode d'emploi
Instrucciones de manejo

Copyright:

**Enraf
Nonius**




Enraf-Nonius B.V.
P.O. Box 12080
3004 GB Rotterdam
The Netherlands
Tel.: +31(0)10 - 20 30 600
Fax: +31(0)10 - 20 30 699




Order number : 1497.750-45
2nd print, August 2005

Sonopuls 490 / 491

Bedieningshandleiding
Operating Instructions
Gebrauchsanweisung
Mode d'emploi
Instrucciones de manejo

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TABLE OF CONTENTS

Table of contents	4
Introduction	19
Terminology	20
Preliminary notes	20
General	20
Installation	20
Connection	20
Technical maintenance	20
Current density	20
Connection and disconnection reactions	21
Electrolytic effects	21
Electromagnetic interference	21
Product Liability	21
International safety standard IEC 601-1	21
The therapy possibilities in a nutshell	21
Ultrasound therapy	21
Two pole interferential current	21
Asymmetrical biphasic pulse current (TENS)	21
Medium-frequency interrupted direct current (Gi)	22
Combination therapy	22
Controls	22
Installation	24
Treatment head as cathode	24
Indications and Contra-indications	24
Operation	25
Operation - Programs	26
Remarks on operation	27
The treatment head: a precision instrument	27
The contact medium	27
Broken contact in the CC-mode	27
Variations in intensity/power	27
Maintenance by the user	27
Environment Information	28
Fault conditions	28
Specifications	29
Technical data	30
Ordering data	30

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General

The Sonopuls 490 and 491 are portable microcomputer-controlled units for ultrasound therapy. The large display, with surrounding the operating keys, adds to a clear display of the parameter settings. The Sonopuls 491 offers both ultrasound therapy and electrotherapy, which can be combined. The units are battery powered and can thus operate independently of the mains power supply. The light weight and compact dimensions of the unit make it suitable for many applications, e.g. for the treatment of patients at home.

We would like to draw your attention to the fact that all text parts regarding to electrotherapy and combination therapy are not applicable for the Sonopuls 490.

Therapeutic possibilities

By using continuous and pulsed ultrasound (1 or 3 MHz as desired) it is possible to achieve muscle relaxation and pain reduction and also to improve the blood perfusion and recovery ability of tissue.

The Sonopuls 491 offers a choice of three current types:

- two-pole interferential current: amplitude modulated medium-frequency alternating current;
- asymmetrical biphasic pulsed current (TENS);
- medium-frequency direct current (duty cycle 95%), strongly galvanic character).

Depending on the selected current type, electrotherapy can be applied for muscle stimulation, pain control by stimulation of nerve tissue and iontophoresis.

As the output characteristic of the electrotherapy channel can be switched, at choice, between Constant Current (CC) or Constant Voltage (CV), the Sonopuls 491 is also suitable for dynamic electrode applications.

The use of the combination key allows any desired combination of ultrasound and electrotherapy current type to be selected.

The therapeutic possibilities include: localisation and treatment of trigger points, hyperaesthetic areas, dermatomes and myotomes.

(For further information, please refer to the therapy books, which can be ordered separately).

Programs

The Sonopuls offers the possibility of storing frequently used settings (maximum 10) in the computer memory. These programs remain stored, even if there is an interruption in the main electricity supply. A special memory (program "0") facility can ensure that, after switching on, the unit automatically sets itself to the most frequently used current type.

Design

In the design of the Sonopuls, more attention than ever has been paid to ease of operation, appearance, functionality and safety.

The result is a unit that combines a logical operating sequence, supported by a clear display, with a modern appearance and compact dimensions. The wide range of possible applications, together with the assurance of optimum safety for patient and user (the Sonopuls fully complies with the IEC 601-1 standards) make the Sonopuls a product of the highest quality.

The battery

The Sonopuls is a battery-powered unit. Ensure that the battery is always adequately charged so that the unit is always ready for use, even in unexpected situations. The working life of the battery is maximised if the unit is always stored with an adequately charged battery. The battery is maintenance free. The Sonopuls was designed to be battery operated during the whole day, under normal circumstances. We advise you to switch off the unit when it is not used. If mains power is available, use the unit as far as possible with the charger connected, as this saves the battery. During charging of the battery, the unit can be used.

Multi-frequency treatment head

The Multi-frequency treatment head of the Sonopuls can supply both 1 and 3 MHz ultrasound. The heads are available with a large and small treatment surface. These treatment heads are also suitable for subaquatic treatments. The excellent beam characteristics, ergonomic design and effective contact control of the multi frequency treatment heads make an optimal treatment possible.

Finally

You have made a wise choice in selecting the Sonopuls 490 or Sonopuls 491. We are confident that your unit will continue to give satisfaction over many years of use. Nevertheless, if you have any queries or suggestions, please contact us or our representative / distributor in your area.

In the book "Electrotherapeutic Terminology in Physical Therapy" (section on clinical Electrophysiology, American Physical Therapy Association, March 1990) a start has been made on the standardisation of Electrotherapeutic terminology. Their recommendations have, as far as possible, been followed in the literature associated with the Sonopuls 491. For the sake of clarity, some of the new terms are given below, together with the (old) terms previously used.

Terms used

- Amplitude : Intensity
- Phase duration : Pulse duration
- Phase interval : Pulse interval
- Direct current : Galvanic current
- Frequency modulation : Spectrum frequency

PRELIMINARY NOTES

General

It is important that you read these operating instructions carefully before using the Sonopuls. Please make sure that these instructions are available to all personnel who operate the equipment.

For the various applications of this apparatus, please refer to our therapy books "Low and Medium Frequency Electrotherapy", reference number 1480.762 and "Ultrasound Therapy", reference number 1482.762.

Pay attention to the following before using the Sonopuls:

1. Keep yourself informed of the contra-indications (see pages 24, 25).
2. The apparatus may not be used in close proximity (i.e. less than 2 metres) to shortwave or microwave equipment.

The manufacturer cannot be held responsible for the results of using this apparatus for any purposes other than described in these operating instructions.

Installation

- Do not install the unit in a location near to a heat source such as a radiator.
- Avoid exposure to direct sunlight, rain, excessive dust, dampness, mechanical vibrations and shocks.
- Cordless telephones should not be used in close proximity of the unit.
- This unit should not be used in so-called "wet-rooms" (hydrotherapy rooms). The unit has to be installed in such way that liquid cannot enter.

Should any liquid ingress into the housing of this equipment, unplug the unit from the wall socket (if connected) and have it checked by an authorised expert.

Connection

- With regard to the use of the sponge pads, please make sure that they do not come into contact with non-intact skin and that they are well moistened with lukewarm water to avoid undesired current densities irritate the skin and causing pain.
- Mains supply connections must comply with the national requirements regarding medical rooms.
- Before connection of this apparatus to the supply, check that the mains voltage and frequency stated on the type plate of the battery charger corresponds to that of the mains supply.
- The battery charger is a part of the supply circuit on which the safety of the equipment is partly based. As a consequence, the approvals of the Sonopuls 490 or Sonopuls 491 are valid in combination with this particular type of battery charger only (ENA-1810).

Do not, therefore, attempt to use any other charger with the Sonopuls 490 or Sonopuls 491.

Technical maintenance

We recommend to have the unit checked annually. This may be done by your supplier, or by another agency, authorised by the manufacturer. It is also recommended that a record of the service history is kept for all activities relating to service and maintenance. In some countries this is even obligatory.

Maintenance and all repairs should only be carried out by an authorised agency. The manufacturer will not be held responsible for the results of maintenance or repairs by unauthorised persons.

Opening of the equipment by unauthorised agencies is not allowed and will terminate any claim to warranty.

Current density

The particular standard for Therapeutic Electrical Stimulation Equipment is IEC 601-2-1 0 (BS5724 part 2.1 0). The maximum permissible Current Density is 2mA R.M.S. per cm². To verify that this value is not exceeded in electrotherapy, divide the current output in mA by the effective electrode area in cm²; a value of less than two should result.

Note that placing the electrodes should be done carefully and solidly, and that the use of small electrodes in combination with high intensities may cause skin irritations or burns.

Connection and disconnection reactions

Constant Current (CC) output characteristics may cause unpleasant connection and disconnection reactions if the electrodes are not securely placed or lose contact with the skin. Make sure the output is set to 0 mA when you apply or remove the electrodes. Preferably use Constant Voltage (CV) in case of dynamic application techniques.

Electrolytic effects

Electrolysis occurs under the electrodes when current types with a DC component are applied. Because the greatest concentration of electrolytic by-products caused by ion migration occur under the electrodes, we recommend the use of the supplied sponges to keep the effects to a minimum. Make sure that the sponges are kept well moistened and place the thick side of the sponge between the flexible rubber electrodes and the patient.

Electromagnetic interference

Simultaneous connection of a patient to High Frequency Therapy (shortwave or microwave) equipment and an electrical stimulator must not be done as this may result in burns at the site of the stimulator electrodes. Therefore, simultaneous connection is not allowed.

Operation in close proximity (i.e. less than 2 metres) to working shortwave or microwave therapy units may produce instability in the Sonopuls output. To prevent electromagnetic interference, we strongly advise using separate mains groups (phases) for the Sonopuls and the shortwave or microwave equipment, and a distance of at least 2 metres is kept between the Sonopuls and the shortwave or microwave equipment. Ensure that the mains cable of the shortwave/microwave unit does not come near the Sonopuls or the patient.

If problems with electromagnetic interference persist, please contact your authorised distributor.

Product Liability

A law on Product Liability has become effective in many countries. This Product Liability law implies, amongst other things, that once a period of 10 years has elapsed after a product has been brought into circulation, the manufacturer can no longer be held responsible for possible shortages of the product.

International safety standard IEC 601-1

This equipment meets the requirements of the International Standard for the safety of medical electrical equipment IEC 601-1 (general standard), IEC 601-2-5 (standard for ultrasound equipment) and IEC 601-2-10 (standard for electrotherapy equipment / only Sonopuls 491).

Ultrasound therapy

Frequency

The Sonopuls offers a choice of two ultrasound frequencies: 1 and 3 MHz. 3 MHz sound waves are absorbed relatively strongly by the tissue and consequently have a slightly less deep effect than 1 MHz waves.

Ultrasoundmode

Both continuous and pulsed ultrasound can be applied. The advantage of the pulsing sound beam is that the thermal sensation is suppressed. In this way it is possible to set a higher intensity than in the case of continuous ultrasound.

Pulse frequency / duty cycle

The pulse frequency is set at 100 Hz. The duty cycle (= average effective power or ratio pulse duration: pulse interval) can be altered in 5, 10, 20, 50, 80 and 100 %.

Two pole interferential current

The Sonopuls 491 uses two pole interference, enabling the interferential current to be generated in the equipment itself, and thus only two electrodes are needed instead of four.

Carrier wave frequency

This parameter is adjustable in 2, 4 or 10 kHz. Depending on the type of nerve fibre (motor/sensory) a low or high carrier wave frequency may be selected.

Stimulation frequency (AMF)

The stimulation frequency can be adjusted in steps between 0 and 150 Hz. The choice of setting is dependent on the nature, stage, severity and localisation of the disorder.

Frequency modulation (spectrum)

This parameter can be adjusted in steps between 0 and 100 Hz, and serves to prevent accommodation.

Frequency modulation programs

There are two frequency modulation (sweep) programs: 6/6 s and 1/30/1/30 s.

Asymmetrical biphasic pulse current (TENS)

Phase duration

The phase duration can be adjusted in steps between 20 and 400 µs. A short phase duration allows selective stimulation of the thick fibres.

Stimulation frequency

The frequency can be adjusted in steps between 5 and 200 Hz.

Burst frequency

Burst TENS is particularly suitable for the treatment of deeper-lying painful regions (myofascial pain) and chronic pain. The burst frequency is adjustable between 1 and 9 Hz.

Frequency modulation (spectrum) and frequency modulation programs

These parameters are also adjustable within TENS. As within the interferential current, accommodation of tissue will be prevented by using frequency modulations and a modulation program. The frequency modulation can be adjusted in steps between 0 and 180 Hz. The next two frequency modulation programs are available: 6/6 s and 1/30/1/30 s.

Medium-frequency interrupted direct current (Gi)

This type of current has a high pulse frequency of 8000 Hz and a duty cycle of approx. 95%. In this way the sensation of the current is experienced less aggressively than continuous direct current, whereas the effects are the same.

Combination therapy

Depending on the depth effect desired and the stage of the disorder, a combination of ultrasound and a current type is selected. Within the Sonopuls 491, it is possible to combine ultrasound with the available current types.

UNIT

(see fold-out sheet at front)

[1] ON/OFF switch

The switch is used to switch the Sonopuls on and off. The unit operates on the battery and/or battery charger.

[2] Type number/warning sticker

Provides information on the unit, such as type and serial number, as well as connection data such as mains voltage and maximum current consumption.

[3] Connector for battery charger

For the connection of the battery charger

[4] Connection patient cable electrotherapy

Connection for the patient cable supplied with the unit. This connection can also be used to connect a Vacotron 560 unit for use with vacuum electrodes.

For combined applications use only Enraf-Nonius type BF equipment. The very low leakage current of this type of equipment ensures absolutely safe therapy.

[5] Connection treatment head

Connecting point for a large or small multi frequency treatment head 1 and 3 MHz.

Connection of treatment heads other than the Enraf-Nonius treatment heads can affect the safety of the patient and correct functioning of the equipment, and is therefore not permitted.

[6] Battery charge indicator

The lamp provides a general indication of the battery charge, there are four indications. The lamp:

- lights continuous green : battery fully charged;
- lights flashing green : battery partial discharged;
- lights yellow * : battery insufficiently charged;
- is out : battery fully discharged, unit is not switched on.

* The current output declines slowly, and the unit switches off automatically. This means that the treatment cannot be continued. The battery must be charged as described in chapter "Maintenance by the user" of this manual.

[7] Charger indicator lamp

The lamp indicates that the unit is connected to the battery charger and mains power supply. The lamp:

- lights green: the battery charger is functioning correctly;
- lights yellow: an incorrect battery charger is connected to the unit.

Connection of a battery charger other than the prescribed charger can adversely affect the safety of the patient and the functioning of the unit, and is therefore not permitted.

CONTROL PANEL

(see fold-out sheet at front and back)

Selection of the therapy form takes place through the push-buttons [8] and [9]. After which the parameter settings are chosen with push-buttons [10] up to and including [14]. The memory is operated through push-button [15] and [16].

The symbol of the chosen setting will blink during a period of 10 sec. after the push-button is pressed. Pressing the push-button again, when the symbol is still blinking, will activate the next possibility as the actual setting chosen (and will blink).

The parameter values of push-button [10] - [14] and [16] are adjusted by means of the Central control [17]. The parameter values can be altered as long as the symbol is blinking.

[8] Selector button therapy form

Using this button, you can select from:

Electrotherapy

Combination therapy

Ultrasound therapy

[9] Selector button electrotherapy current type

A selection can be made of the following current types:

Two-pole interferential current

Asymmetrical biphasic pulse current (TENS)

Medium-frequency interrupted direct current (Gi)

[10] Selector button carrier wave frequency, phase duration and ultrasound frequency

Depending on the chosen therapy form, the following parameter settings can be changed:

- carrier wave frequency (kHz) with two-pole interferential current
- phase duration (μ s) with asymmetrical biphasic pulse current
- ultrasound frequency (MHz) with combination therapy and ultrasound therapy

[11] Selector button stimulation frequency and duty cycle ultrasound

Button for the adjustment of the stimulation frequency for two-pole interferential current, asymmetrical biphasic pulse current and duty cycle for ultrasound therapy.

[12] Selector button memory - Sonopuls 490

By means of this push-button a choice can be made to store or retrieve a program from the memory. The program number is changed by the Central control [17]. For the operation of the memory, we refer to chapter "operation" of this manual.

[12] Selector button frequency modulation and modulations programs - Sonopuls 491

Button for the adjustment of the frequency modulation (spectrum). As soon as a frequency is set the modulation-symbol will become visible on the display. After pressing the push-button once more, the modulation-symbol will blink and a modulation program can be chosen with the Central control [17].

[13] Selector button burst frequency and timer

Adjustment of a burst frequency is only possible for asymmetrical biphasic pulse current (and when no frequency modulation is set).

[14] Selector button CC/CV and display function ultrasound

The electrotherapy channel can be used with either constant current (CC) or constant voltage (CV). The CV characteristic is recommended for use with dynamic electrode techniques, and for combination of electrotherapy with ultrasound.

The following parameters can be indicated on the display:

- ultrasound intensity in W/cm^2
- average power in Watt

[15] Confirm button

Button for the confirmation of "to retrieve a program from" or "store a program into" the memory.

[16] Memory button - Sonopuls 491

With this push-button a choice can be made to retrieve a program from the memory or to store a program into the memory. The program number is changed by the Central control [17]. For the operation of the memory, we refer to chapter "operation" of this manual.

[17] Central control

Regulator for adjusting the parameter values.

[18] Amplitude control

Regulator without end-stop for setting the patient current.

[19] Intensity control

Regulator without end-stop for setting the ultrasound intensity.

MULTI-FREQUENCY TREATMENT HEAD**[20] Contact control indicator**

The indication ring of the selected treatment head lights up. In case of sufficient acoustic contact, the ring extinguishes and ultrasound power is generated.

INSTALLATION

Connection of the battery charger

- Connect the supplied battery charger to the connector for battery charger [3].
- Connect the battery charger to a wall socket.

Charging before the first use of the unit

The Sonopuls is delivered with a fully charged battery. We nevertheless advice to charge the battery for a period of at least 12 hours before the Sonopuls is used for the first time.

Switching on and self test

- Switch on the unit using the ON/OFF switch [1].
- Immediately after switching on, the unit carries out a self test.

Check that:

- All LED's light up briefly.
- At the end of the test a beep can be heard.

Connecting the accessories

- Connect the supplied patient cable to connection [4].
- Connect the supplied multi-frequency treatment head to connection [5].

TREATMENT HEAD AS CATHODE

The treatment surface of a connected treatment head is always connected with the negative electrode (cathode) of the electrotherapy channel.

Combination therapy:

- Remove the negative electrode (black plug) of the electrode cable.
- The current flows between the positive electrode (red plug) and the treatment head.

Electrotherapy:

- Under no circumstances should the treatment surface of the ultrasound treatment head come into contact with the electrodes in the patient circuit whilst the unit is operating.

INDICATIONS AND CONTRA-INDICATIONS

Please refer to the optional therapy books.

Indications (Ultrasound)

- conditions affecting joints, bone, muscle tissue;
- rheumatoid arthritis (not acute);
- affections of the peripheral nerves;
- affections of the blood circulation;
- affections of the internal organs;
- affections of the skin, scar tissue;
- Dupuytren's contracture;
- open wounds, decubitus ulcer, post-traumatic conditions.

Absolute specific contra-indications (Ultrasound)

- eyes;
- heart;
- pregnancy;
- epiphyseal discs;
- brain tissue;
- testicles.

Relative specific contra-indications (Ultrasound)

- status post laminectomy;
- loss of sensation;
- endoprostheses;
- tumours;
- post-traumatic sequelae;
- thrombophlebitis and varices;
- septic inflammation;
- diabetes mellitus.

Indications (electrotherapy)

Diagnostics, electro palpation for:

- pain points;
- trigger points;
- hyperaesthetic regions;
- motor trigger points;

Therapy:**A. Pain control in..**

- pain points;
- trigger points;
- hyperaesthetic regions.

B Autonomic disorders such as:

- shoulder/hand syndrome;
- Raynaud's disease;
- Buerger's disease;
- Sudeck's dystrophy;

C. Muscle stimulation:

- for restoration of kinaesthesia;
- for treatment of atrophy;
- for stimulation of the internal and external sphincters in the treatment of urinary incontinence;
- revalidation;
- muscle strengthening (sport);

D. Iontophoresis.**Contra-indications (electrotherapy)**

- pyrexia;
- tumours;
- tuberculosis;
- localised inflammation;
- thrombosis;
- pregnancy;
- pacemakers;
- metal implants.

Indications (combined therapy)

- myofascial trigger points;
- referred sensation;
- incipient neurogenous inflammation.

Contra-indications (combined therapy)

- As for ultrasound and electrotherapy.

OPERATION

General

Switch the unit on [1]. The display shows the program settings as programmed under program number 0. A preferred program has been entered in the factory. This can be replaced by your own preferred settings. For the procedure, we refer to paragraph "Programs" of this chapter.

Select with push-button [8] the desired therapy type. Blinking of the symbol indicates the chosen therapy type. The adjustment will change by operating the push-button once again.

Select with push-button [9] the desired current type (if applicable).

The different parameters are selected with the push-buttons [10] / [14]. One push-button operates all parameters which lies on the same horizontal line. Parameters are set by means of the central control [17]. As long as the symbol of the parameter is blinking the parameter can be altered. The display only displays the parameters which are applicable for the chosen therapy form.

Operation Sonopuls 490:

- Connect the ultrasound treatment head.
- Operate button [10] once, the "frequency" symbol and MHz-indication blink [I], set the ultrasound frequency by means of the central control [17].
- Operate button [11] once, the "duty cycle" symbol blink [K], adjust the desired duty cycle by means of the central control [17].
- Operate button [13] once, the "min" symbol blink [Q], set the treatment time by means of the central control [17].
- Operate button [14] once, the "display mode" symbol blink [S], alter the adjustment by means of the central control [17].
- Set the ultrasound intensity by means of the intensity control [19].
- Position the ultrasound treatment head on the patient (using gel as a contact medium).

Operation Sonopuls 491:

- Connect the patient cable and ultrasound treatment head.
- Use only the positive electrode of the patient cable. This is identified by a red connector sleeve.
- Select with button [8] "combination therapy" [B].
- Select with button [9] "two-pole interferential therapy" [D].
- Operate button [10] once, the "frequency" symbol and kHz-indicator [G] blink, alter the carrier frequency by means of the central control [17].
- Operate button [10] once, the "frequency" symbol and MHz-indicator [I] blink, set the ultrasound frequency by means of the central control [17].
- Operate button [11] once, the "stimulation frequency" symbol [J] blink, set the frequency by means of the central control [17].
- Operate button [11] once, the "duty cycle" symbol [K] blink, set the duty cycle by means of the central control [17].
- Operate button [12] once, the "frequency modulation" symbol [L] blink, adjust the frequency by means of the central control [17]. As soon as the frequency is set the "frequency modulation programs" symbol [M] appears on the display.
- Operate button [12] once, the "frequency modulation programs" symbol [M] blink, alter the modulation program by means of the central control [17].
- Operate button [13] once, the "min" symbol [Q] blink, set the treatment time by means of the central control [17].

- Operate button [14] once, the “output characteristic” symbol [R] blink, alter the adjustment by means of the central control [17].
- Operate button [14] once, the “display mode” symbol [S] blink, alter the adjustment by means of the central control [17].
- Set the ultrasound intensity by means of the intensity control [19].
- Position the ultrasound treatment head on the patient (using gel as a contact medium).
- As soon as adequate acoustic contact has been established, the current amplitude can be increased [18].

OPERATION - PROGRAMS

General

Frequently used equipment settings can be stored under the program numbers 0 to 9. The program numbers are shown on the display, after operation of button [16] (Sonopuls 491) or [12] (Sonopuls 490).

With the exception of the current amplitude, all setting can be programmed, including the frequencies of the ultrasound treatment heads connected.

Entering a program in the memory - Sonopuls 490

- Set the required settings which need to be programmed.
- Operate button [12] once, the “memory retrieve” symbol [N] blink, select, by means of the central control [17], a program number (the display will show the settings which are now stored under that number).
- Operate button [12] once, the “memory store” symbol [O] blink. The display will now show your new settings which need to be stored.
- Operate button [15] to confirm the storage of the new settings (and overwrite the old settings). The confirmation must be done before the “memory store” symbol stops to blink.

Entering a program in the memory - Sonopuls 491

- Set the required settings which need to be programmed.
- Operate button [16] once, the “memory retrieve” symbol [N] blink, select, by means of the central control [17], a program number (the display will show the settings which are now stored under that number).
- Operate button [16] once, the “memory store” symbol [O] blink. The display will now show your new settings which need to be stored.
- Operate button [15] to confirm the storage of the new settings (and overwrite the old settings). The confirmation must be done before the “memory store” symbol stops to blink.

Retrieving a program from the memory - Sonopuls 490

- Operate button [12] once, the “memory retrieve” symbol [N] blink, select, by means of the central control [17], a program number.
- Operate button [15] once to confirm your choice.

Retrieving a program from the memory - Sonopuls 491

- Operate button [16] once, the “memory retrieve” symbol [N] blink, select, by means of the central control [17], a program number.
- Operate button [15] once to confirm your choice.

Interrupting the procedure

If the program button have been operated (deliberately or accidentally) the procedure can be interrupted by pressing any other button, with the exception of the program button or confirm button.

Program 0; a special program

After the unit has been switch on, and the self-test has been completed, program 0 is selected automatically. A preferred program has been entered in the factory. This can be replaced by your own preferred settings by following the above procedure. The unit will then automatically start by selecting your preferred program.

The treatment head: a precision instrument

A treatment head is a precision instrument. Great care is taken in the development and production in order to obtain the best possible beam characteristics. Rough treatment (jarring or dropping) can adversely affect these characteristics, and must therefore be avoided.

The contact medium

In order to ensure efficient transfer of energy, a contact medium is required between the treatment head and the body. Air causes virtually total reflection of the ultrasound energy. The best medium for the transfer of ultrasound energy is a gel.

- For preference, use Enraf-Nonius Contact-Gel®, as this allows the excellent characteristics of the treatment heads to be used to their full advantage.
- The gel should be applied to the part of the body to be treated and then spread out with the treatment head.
- Never apply the gel to the treatment head. The treatment head will register this as contact and may emit ultrasound energy, which could damage the treatment head.

If the body surface is very irregular, making it difficult to obtain good contact between the treatment head and the body, or if direct contact must be avoided (e.g. due to pain), the affected area may be treated under water (subaqual method). The water should be degassed (by previous boiling) in order to prevent air bubbles arising on the treatment head and the body.

Enraf-Nonius Contact-Gel® is a registered trademark

Broken contact in the CC-mode

When contact is broken in the CC-mode, the equipment will indicate this with a beep-tone. In order to avoid an unpleasant sensation for the patient, when the treatment head or electrode is placed back onto the surface that is to be treated, the current will ramp down to 0. To continue the treatment, you will have to set the current again.

Variations in intensity/power

The Sonopuls contains a control circuit which ensures that the ultrasound intensity in the tissue (= wave motion of the "micro-massage") is kept at a constant level as far as possible.

In three situations there is a deviation in intensity as indicated on the display of the Sonopuls:

Situation 1: Incomplete contact

When the treatment head only contacts partially with the body, the set intensity is automatically reduced to prevent too much energy being introduced into the patient's body via the small surface. When less than 80% of the surface of the treatment head makes

contact with the patient the intensity is interrupted (=contact control). The display starts flashing and the treatment time is interrupted.

Situation 2: Differences in tissue density

The Enraf-Nonius treatment heads are factory-set to issue energy in water (according to the IEC- standards). However, body tissue reacts differently to ultrasound than water; it is harder to set body tissue into vibration. Furthermore, the various tissue structures of the body (fat tissue, muscles, bone structures) all react differently to ultrasound. In case of a certain intensity setting the Sonopuls automatically adapts the supplied energy to maintain the wave motion of the micro-massage (in the various tissues). During this treatment the display indicates this intensity variation.

Situation 3: Reflections

Bone structures in particular reflect ultrasound. These reflections result in an increased intensity and an intensification of the micro-massage. The increased intensity is limited at a maximum of 20% of the set value.

A maximum safety limit ensures that the intensity can never exceed 3 W/cm² specified in the IEC standard.

MAINTENANCE BY THE USER

Apparatus

Before performing any user maintenance, switch off the apparatus and disconnect it from the power supply. The apparatus can be cleaned with a damp cloth. Use lukewarm water and a non-abrasive liquid household cleaner (no abrasive, no alcohol content solution).

Charging

We advise charging the battery overnight. The charging time will then generally be more than adequate. However, in the case of a completely exhausted battery (lamp [6] is out), the charging time can be longer, e.g. 15 - 20 hours. The Sonopuls can be used during the charging period. Switching off of the Sonopuls will shorten the charging time. Ensure that the charger indicator lamp [7] is lit, indicating that the charger is functioning correctly.

Treatment heads

The treatment heads and cables should be regularly inspected for damage, e.g. hairline cracks, which could allow penetration by liquids. Clean the contact surface immediately after each treatment. Make sure that no ultrasound gel remains on the treatment head. We further recommend cleaning the head and cable daily, using lukewarm water. The treatment heads can be disinfected using a cloth moistened with 70% alcohol.

Electrodes and sponge pads

The rubber electrodes should be cleaned with luke-warm water. For stubborn stains or dirt, cleaning and disinfection of the rubber electrodes can be done using a 70% alcohol solution. During this process, the black colouration of the rubber electrodes may be stained. This does not affect the operation of the electrodes.

The sponge pads must be washed after use in warm water with a household cleaner. Then they must be rinsed with water, drained thoroughly and then dried, after they have been placed like roofing tiles. Damaged sponge pads must be replaced.

In areas where the water is soft, it is possible that the water may conduct poorly resulting in abnormally low currents. In this case it is possible to use a saline solution to improve the conductivity of the water.

When the electrodes are not being used, the sponge pads should be removed. This increases the operational life of the rubber electrodes. We advise that you should keep an extra set of sponge pads in stock.

Patient cable

Check the cable regularly for damages and/or bad electrical contact. We also advise that you keep an extra patient cable in stock.

Environment Information

Your Sonopuls contains materials which can be recycled and/or are noxious for the environment. At the end of the duration of life, specialised concerns can take apart the Sonopuls and sort out the noxious materials and materials for recycling. By doing so you contribute to a better environment.

Please ensure that you are well informed of the local rules and regulations regarding to the removal of equipment and accessories.

Lamps and display fail to light up

- The battery is possibly exhausted. Try to use the Sonopuls with the battery charger connected to the mains supply.

Error message “Error code ...” on the display

The unit has discovered a fault using its “self-test”. Remove any applicators or cables from the output sockets from the equipment, switch the unit back on. If the code re-appears, contact your supplier. The unit is probably defective.

The contact control fails to operate

Check that the treatment head has been selected. Check the contact control by placing the treatment head in a bowl of water. It may be that the contact control is working, but that the contact gel is inadequate. If several treatment heads are available, it is possible to determine whether the problem lies in the gel, the treatment head or the unit itself.

Insufficient or no output current

- Check the patient cable for breaks or poor contact.
- Ensure that the sponges are sufficiently damp.
- If necessary, use a salt solution.
- Clean the electrodes and check for poor contact.
- In combination therapy, the maximum output current for two-pole interferential therapy is limited to 60 mA.

Ultrasound

Frequencies	: 1 and 3.1 MHz
Display	: intensity in W/cm ² and power output in W (SATP*)
Contact control threshold	: 80%
Timer	: 0 - 30 minutes, coupled to contact control

Ultrasound, continuous

Pulse frequency / duty cycle	: 100 Hz / 100 %
Number of connections	: 1
Intensity	: 0 - 1.5 W/cm ²

Ultrasound, pulsed

Pulse frequency / duty cycle	: 100 Hz / 5, 10, 20, 50, 80 %
Number of connections : 1	
Intensity	: 0 - 3 W/cm ² , duty cycle 5, 10, 20, 50 % 0 - 1.8 W/cm ² , duty cycle 80 %

Treatment heads

1 and 3 MHz, large (standard)	: Geometric surface 5.8 cm ² ERA* 5 cm ² , BNR* max. 5 type collimating, side panel radiation max. 10 W/cm ² (IEC requirement < 100 mW/cm ²)
1 and 3 MHz, small (option)	: Geometric surface 1.15 cm ² , ERA* 0.8 cm ² , BNR* max. 5 type collimating, side panel radiation max. 10 mW/cm ² (IEC requirement < 100 mW/cm ²)

*SATP = Spatial Average Temporal Peak (average pulse power)

*ERA = Effective Radiating Area, this is the effective radiating area of the treatment head.

*BNR = Beam Nonuniformity Ratio, indicates the ratio between the peaks and the average value of the intensity in the sound beam. A low BNR excludes high, unwanted energy concentrations.

ERA and **BNR** have been measured according the FDA-method. All treatments heads have been provided with contact control.

Current general

Number of channels	: 1
Output characteristics	: Constant Current (CC) or Constant Voltage (CV)
Resolution of CV or CC display	: in steps of 1 mA
Timer	: 0 - 30 min.

The maximum amplitude is achieved within the specification at a load of max. 500 Ohm (CC).

Two-pole interferential therapy

Carrier frequency	: 2, 4 or 10 kHz
Stimulation frequency (AMF)	: 0-5-10-15-20-30-40-50-60-80-100-150 Hz
Frequency modulation (spectrum)	: 0-10-20-30-40-50-80-100 Hz
Frequency modulation programs	: 6/6s and 1/30/1/30s
Amplitude	: 0 - 100 mA 0 - 60 mA in combination therapy

Asymmetrical biphasic pulsed current

Phase duration	: 20-30-40-50-75-100-125-150-175-200-300-400 µs
Frequency	: 5-10-15-20-30-40-50-60-80-100-120-140-160-180-200 Hz
Burst frequency	: 1, 2, 3, 4, 5, 6, 7, 8, 9 Hz
Frequency modulation (spectrum)	: 0-5-10-20-30-40-50-60-80-100-120-140-160-180 Hz
Frequency modulation programs	: 6/6s and 1/30/1/30s
Amplitude	: 0 -140 mA, if phase duration ≤ 200 µs 0 - 100 mA, if phase duration is 300 µs 0 - 70 mA, if phase duration is 400 µs

Medium-frequency direct current

Frequency	: 8000 Hz
Duty cycle	: 95%
Amplitude	: 0-40 mA

Programming

Memory	: memory for 10 sets of freely selectable therapy parameters
--------	--

Battery charger

Model number	: ENA-1810
Mains voltage	: 230 V EUR 230 V BSI (optional) 120 V UL/CSA (optional)
Frequency	: 50/60 Hz
Permissible variation	: $\pm 10\%$
Current consumption	: 17 VA
Output voltage	: 18 V
Max. Output current	: 500 mA

Sonopuls

Battery	: 12 V, 1.8Ah maintenance-free lead accumulator; dim.: 178.5 x 34 x 60.5 mm (l x w x h)
Patient leakage current	: typically 1 μ A (requirement IEC $\leq 100 \mu$ A)
Ditto, first fault condition	: typically 2 μ A (requirement IEC $\leq 500 \mu$ A)
Safety class	: II* type BF**, according to IEC 601-1
Weight	: 4.0 kg
Dimensions	: 290 x 280 x 105 mm
Safety tests	: TÜV Rheinland

Environment conditions for Transport and Storage

Environment temperature	: -10° till +50° C
Relative humidity	: 10 till 100 %
Atmospheric pressure	: 500 till 1060 hPa

This equipment complies with all requirements of the Medical Device Directive (93/42/EEC), Class IIb.

*II : Safety class II (double insulated).

**BF : The equipment has a floating patient circuit. The Sonopuls 490 and Sonopuls 491 meets the requirements of IEC 601 -1, IEC 601-2-5 and IEC 601-2-1 0.

Copies of test reports available on application.

Technical modifications reserved.

ORDERING DATA

For the ordering data of the Sonopuls, standard accessories and additional accessories we refer to the Catalogue for Physiotherapy.

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