
Sonopuls 590

Bedieningshandleiding
Operating Instructions
Gebrauchsanweisung
Mode d'emploi
Instrucciones de manejo



Copyright:

**Enraf
Nonius**

Enraf-Nonius B.V.
P.O. Box 12080
3004 GB Rotterdam
The Netherlands
Tel.: +31(0)10 - 20 30 600
Fax: +31(0)10 - 20 30 699

Article number 1495.753-42
2nd print, August 2005

Sonopuls 590

Bedieningshandleiding
Operating Instructions
Gebrauchsanweisung
Mode d'emploi
Instrucciones de manejo

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TABLE OF CONTENTS

Illustration of Sonopuls 590	2
Foreword	12
Preliminary notes	12
General	12
Installation	12
Connection	12
Electromagnetic interference	12
Technical maintenance	11
Product liability	11
International safety standard IEC 601-1	11
Radio interference suppression	11
Indications and contra-indications	11
Controls	12
Operating sequence	13
Remarks on operation	14
The treatment head; a precision instrument	14
The contact medium	14
Variations in intensity	14
Maintenance by the user	14
Fault conditions	14
Specifications	15
Technical data	15
Ordering data	15

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General

The Sonopuls 590 is a microprocessor controlled unit for continuous and pulsed ultrasound therapy. The equipment is only meant to be used by competent personnel in physiotherapy, rehabilitation or adjacent areas. A great deal of attention has been paid to design, functionality and safety. The carefully chosen dimensions and the built-in carrying handle also make the unit suitable for ambulant use.

Two ultrasound frequencies

Therapeutic effects are only obtained when the applied energy is absorbed in the tissue. For deeper-lying tissues this occurs at a frequency of 1 MHz, while 3 MHz is very suitable for superficial structures. Large and small treatment heads are available for both frequencies.

Two treatment head connections

The Sonopuls 590 allows two treatment heads to be connected simultaneously, e.g. a large 1 MHz head and a small 3 MHz head. The treatment head to be used can then be selected from the unit. The treatment heads are fully interchangeable, so that a treatment head can be connected to any other Sonopuls 590 unit without recalibration.

Beam characteristics

The quality of the therapy depends, to a great extent, on the quality of the ultrasound unit and the treatment heads. The characteristics of the ultrasound beam can only be measured with extremely expensive test equipment. For this reason, safety authorities require that characteristics such as the ERA (Effective Radiating Area) and the BNR (Beam Nonuniformity Ratio) should be stated on the treatment head. Random tests are regularly carried out by an independent institute to check the accuracy of the figures given.

Efficiency and safety

All treatment heads are provided with contact control. This feature, in combination with the excellent beam characteristics, maintains the quality of the treatment. Maximum safety is inherent in the design of the unit. Both the electrical safety and the ultrasound safety of the Sonopuls 590 comply fully with the IEC 60601-1 safety standards. A built-in microprocessor monitors correct functioning of the unit on switching on, and during operation.

Finally

You have made a wise choice in selecting the Sonopuls 590. We are confident that your unit will continue to give satisfaction over many years of use. Nevertheless, if you have any queries or suggestions, please contact your authorized distributor.

General

It is important that you read these operating instructions carefully before using the Sonopuls 590. Please make sure that these instructions are available to all personnel who operate the equipment.

For the various applications of this apparatus, please refer to our therapy book "Ultrasound Therapy", reference number 1482.762.

Pay attention to the following before using the Sonopuls 590:

1. Keep yourself informed of the contra-indications (see page 11).
2. The apparatus may not be used in close proximity (i.e. less than 2 metres) to shortwave or microwave equipment.
3. The patient must be in sight of the therapist at all times.

The manufacturer cannot be held responsible for the results of using this apparatus for any purposes other than described in these operating instructions.

Installation

- Do not install the unit in a location near a heat source such as a radiator.
- Avoid exposure to direct sunlight, excessive dust, dampness, mechanical vibrations and shocks.
- This unit should not be used in so-called "wet rooms" (hydrotherapy rooms). This unit has to be installed in such way that liquid cannot enter.

Should any liquid ingress into the housing of this equipment, unplug the unit from the wall socket and have it checked by an authorized expert.

Connection

- Mains supply connections must comply with the national requirements regarding medical rooms.
- Before connection of this apparatus to the supply, check that the mains voltage and frequency stated on the type plate of this apparatus correspond to that of the mains supply. 220-240 V model for countries with mains voltage of 220, 230 or 240 volt. 115 V model for countries with mains voltage of 110, 120 or 130 volt.
- This apparatus may only be connected to an earthed wall socket.
- Only use the mains cable supplied with the apparatus.
- A potential equalisation cable may be connected to this unit, if prescribed in the national regulations.

Electromagnetic Interference

Operation in close proximity (i.e. less than 2 metres) to working shortwave or microwave therapy unit may produce instability in the Sonopuls 590 output.

To prevent Electromagnetic Interference, we strongly advise using separate mains groups (phases) for the Sonopuls 590 and the shortwave or microwave equipment. Ensure that the mains cable of the shortwave/microwave unit does not come near the Sonopuls 590 or the patient.

If problems with Electromagnetic Interference persist, please contact your authorized distributor.

Technical maintenance

We recommend to have the unit checked annually. This may be done by your supplier or by an other agency, authorized by the manufacturer. It is also recommended that a record of the service history is kept for all activities relating to service and maintenance. In some countries this is even obligatory.

The check up and/or technical maintenance must be carried out conform the procedure described in the service manual of the unit.

No attempt should be made to open the unit. Maintenance and all repairs should be carried out by an authorised agency.

The manufacturer will not be held responsible for the results of maintenance or repairs by unauthorized persons.

Product liability

A law on Product Liability has become effective in many countries. This Product Liability law implies, amongst other things, that once a period of 10 years after a product has been brought into circulation has elapsed, the manufacturer can no longer be held responsible for possible shortages of the product.

International safety standard IEC 60601-1

This equipment meets the requirements of the International Standard for the safety of medical electrical equipment IEC 60601 -1, (general standard), IEC 60601-2-5 (standard for ultrasound equipment)

USA only:

Safety Precautions

This unit operates with high voltages. Servicing of the SONOPULS 590 should be referred to qualified service technicians or returned directly to the distributor. To assure continued compliance with FDA 21 C.F.R. 1050.10 standards, the SONOPULS 590 should be calibrated and safety tested at least once each year. The producer recommends that all service work be performed by the distributor. However, a service manual can be obtained by writing to the factory for those who are qualified to service the device and are fully aware of their responsibilities.

Caution:

Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous exposure to ultrasonic energy.

INDICATIONS AND CONTRA-INDICATIONS

Please refer to therapy book supplied.

Indications

- Conditions affecting bone tissue, joints, muscle tissue;
- Rheumatoid arthritis (not acute);
- Affections of the peripheral nerves;
- Affections of the blood circulation;
- Affections of the internal organs;
- Affections of the skin, scar tissue;
- Dupuytren's contracture;
- Open wounds, decubitus ulcer, post-traumatic conditions.

Absolute specific contra-indications

- Eyes;
- Heart;
- Gravid uterus;
- Epiphysial discs;
- Brain tissue;
- Testicles.

Relative specific contra-indications

- Post laminectomy;
- Loss of sensation;
- Endoprostheses;
- Tumours;
- Post-traumatic sequelae;
- Thrombophlebitis and varices;
- Septic inflammation;
- Diabetes mellitus.

APPARATUS

(see fold-out sheet at front)

[1] Type number/warning sticker

Provides information on the unit, such as type and serial number, as well as connection data such as mains voltage and current consumption.

[2] Potential equalization point

For use in areas where local regulations require connection of a potential equalization cable. An earth (ground) cable and plug can be supplied.

[3] Mains switch, mains input and fuses

The switch is used to switch the Sonopuls 590 on and off. The mains cable is connected to the mains input, which is provided with two fuses.

[4] Carrying handle (retractable)**CONTROLPANEL****[5] 100 Hz pulsed mode selector**

Pulsed ultrasound enables a high intensity to be applied while the mean power remains low. The ratio of pulse duration to pulse interval is 1:4 (duty cycle = 20 %).

[6] Continuous mode selector

For continuous ultrasound with a high mean power.

[7] 48 Hz pulsed mode selector

The pulse frequency of 48 Hz is a multiple of 16 Hz. The latter is a fundamental frequency of the calcium system. The ratio of pulse duration to pulse interval is 1:4 (duty cycle = 20 %).

[8] 16 Hz pulse mode selector

This mode combines a low pulse frequency with a high-energy pulse.

16 Hz is a fundamental frequency of the calcium system. The ratio of pulse duration to pulse interval is 1:4 (duty cycle = 20 %).

[9] Treatment time controls

These controls are used to set the treatment time in steps of 1 minute. Under 10 minutes the treatment time can be set in intervals of half a minute using the minus-control.

[10] Treatment head selector

Selector for treatment head A or B.

[11] Intensity controls

These controls are used to set the intensity. The maximum intensity can be set from zero by one touch on the minus-control.

[12] Display function selector (W/cm² or W)

Switches the display function between intensity (W/cm²) and power (W). When the display is in power (W) the unit allows for the area (ERA) of the treatment head.

[13] Treatment head connection A

Suitable for connection of a large or small 1 MHz or 3 MHz treatment head.

[14] 1 MHz indicator lamp

Indicates that a 1 MHz treatment head is in use.

[15] 3 MHz indicator lamp

Indicates that a 3 MHz treatment head is in use.

[16] Treatment head connection B

Suitable for connection of a large or small 1 MHz or 3 MHz treatment head.

Connection of treatment heads other than those supplied for the Sonopuls 590 can adversely affect the safety of the patient and correct functioning of the equipment, and is therefore not permitted.

ILLUMINATED DISPLAYS**[17] Treatment time display**

When there is adequate ultrasound contact, the timer will count down to zero. If the ultrasound contact becomes inadequate, the timer stops automatically.

[18] Intensity display

Depending on the choice of display function, this display shows the intensity (W/cm² or the power (W). The power depends on the size of the treatment head. The setting flashes as long as there is no adequate ultrasound contact.

TREATMENT HEAD**[19] Contact control indication**

When a treatment head is selected, LED's on the treatment head light up. When the ultrasound contact is adequate, the LED's go out, and ultrasound power is applied.

Ultrasound mode	display function	display indicates
Continuous wave	'W'	The average power in the continuous mode
Continuous wave	'W/cm ² '	The average effective intensity in the continuous mode
Pulsed wave	'W'	The peak pulse ultrasonic power in the pulse-modulated mode
Pulsed wave	'W/cm ² '	The peak pulse effective intensity in the pulse-modulated mode

Mains switch, mains input


- Connect the mains cable supplied to the input.
- Connect mains cable to an earthed (grounded) wall socket.

Immediately after switching on, the unit carries out a self-test.

Treatment head connections A and B

- Connect the treatment heads to connections A and B.

Choice of therapy mode

- Select continuous or pulsed ultrasound.

The selection is indicated on the LED.

Treatment time


- Select the treatment time.

Intensity

- Select treatment head A or B.
- Select display in W/cm² or W.
- Set intensity/power.


The display will flash until adequate contact is made.

Starting treatment

- Position the treatment head (do not forget to apply gel).
- When adequate contact is made, the timer will start.

The contact indicator LED's on the treatment head are no longer lit. The display shows the actual intensity or power.

Ending treatment

Treatment can be stopped by removing the treatment head from the patient and setting the timer to zero. Treatment stops automatically when the treatment time has elapsed.

REMARKS ON OPERATION

The treatment head: a precision instrument

A treatment head is a precision instrument. Great care is taken in the development and production in order to obtain the best possible beam characteristics. Rough treatment (jarring or dropping) can adversely affect these characteristics, and must therefore be avoided.

The contact medium

In order to ensure efficient transfer of energy, a contact medium is required between the treatment head and the body. Air causes virtually total reflection of the ultrasound energy. The best medium for the transfer of ultrasound energy is a gel.

- For preference, use Enraf-Nonius Contactgel as this allows the excellent characteristics of the treatment heads to be used to their full advantage.
- The gel should be applied to the part of the body to be treated, and then spread out with the treatment head.
- Never apply the gel to the treatment head. The treatment head will register this as contact, and may emit ultrasound energy.

If the body surface is very irregular, making it difficult to obtain good contact between the treatment head and the body, or if direct contact must be avoided (e.g. due to pain), the affected area may be treated under water (subaquatic method).

The water should be degassed (by previous boiling) in order to prevent air bubbles arising on the treatment head and the body.

Variations in intensity/power

Any variations between the ultrasound energy actually applied (and displayed) with respect to the intensity/power value set, e.g. due to reflections, are automatically kept within acceptable limits. These limits have been chosen such that the variation can never exceed 20 %

A maximum safety limit ensures that the intensity can never exceed 3 W/cm² specified in the IEC standard.

USA only:

Ultrasound emission can be terminated in the following ways:

- Lifting the treatment head from the patient. Because of the loss of acoustic contact the SONOPULS 590 stops emission.
- Decrease the treatment time to zero with the time setting switch.
- Switch the power off by use of the power switch.

MAINTENANCE BY USER

Apparatus

Before performing any user maintenance, switch off the apparatus and disconnect it from the power supply. The apparatus can be cleaned with a damp cloth. Use lukewarm water and a non-abrasive liquid household cleaner.

Treatment heads

Clean the contact surface immediately after each treatment. Make sure that no ultrasound gel remains on the treatment head. We further recommend cleaning the head and cable daily, using lukewarm water.

The treatment heads can be disinfected using a cloth moistened with 70% alcohol.

Environment Information

Your Sonopuls contains materials which can be recycled and/or are noxious for the environment.

At the end of the duration of life, specialised concerns can take apart the Sonopuls and sort out the noxious materials and materials for recycling. By doing so you contribute to a better environment.

Please ensure that you are well informed of the local rules and regulations regarding to the removal of equipment and accessories.

FAULT CONDITIONS

Lamps and displays fail to light up

- Check the mains cable and mains fuses.

Replacement of fuses

There are two mains fuses contained within the mains input connector. Remove the power cord. Access to the fuse holder is now possible. To open the fuse holder, unscrew the fuse holder cap (counterclockwise). Replace the fuse with one of the same rating (1.6A/ type H). Tighten the fuse holder cap back into the fuse holder (clockwise) and reconnect the power cord.

A number (code) appears on the display

The unit has discovered a fault using its "self-test". Remove any applicators, treatment heads or cables from the output or remote control sockets from the equipment, switch the unit off for 5 seconds then switch it back on again. If the code re-appears, contact your supplier. The unit is probably defective.

A number only appears on the display when the unit is switched on with a treatment head connected

The treatment head is probably defective. If other treatment heads are available, this can be checked by exchanging the treatment head.

The contact control fails to operate

Check that the treatment head has been selected. Check the contact control by placing the treatment head in a bowl of water. It may be that the contact control is working, but that the contact gel is inadequate. If several treatment heads are available, it is possible to determine whether the problem lies in the gel, the treatment head or the unit itself.

Treatment unit

Frequencies	:	1 MHz ($\pm 0.2\%$) and 3 MHz ($\pm 0.2\%$)
Ultrasound mode	:	continuous pulsed at 100Hz, 2 ms pulse duration, duty cycle 20 % pulsed at 48Hz, 4 ms pulse duration, duty cycle 20 % pulsed at 16Hz, 12 ms pulse duration, duty cycle 20 %
Intensity	:	0 - 2 W/cm ² continuous and 0 - 3 W/cm ² pulsed
Display	:	intensity in W/cm ² or power in W (SATP*)
Indication accuracy	:	$\pm 20\%$ (for any level above 10 % of maximum)
Timer	:	0 - 30 minutes, coupled to contact control
Contact control threshold	:	80 %

Treatment heads

1 MHz, large (standard)	:	geometrical area 6.2 cm ² , ERA** 5.0 cm ² , BNR*** max. 6 type collimating, parasitic side radiation max. 10 mW/cm ²
1 MHz, small (optional)	:	geometrical area 1.4 cm ² , ERA 0.8 cm ² , BNR max. 6 type diverging, parasitic side radiation max. 10 mW/cm ²
3 MHz, large (optional)	:	geometrical area 6.2 cm ² , ERA 5.0 cm ² , BNR max. 6 type collimating, parasitic side radiation zero
3 MHz, small (optional)	:	geometrical area 0.7 cm ² , ERA 0.5 cm ² , BNR max. 6 type diverging, parasitic side radiation zero

* SATP= Spatial Average Temporal Peak (mean pulse power).

** ERA = Effective Radiating Area, i.e. the effective area from which the ultrasound beam is emitted.

***BNR = Beam Nonuniformity Ratio, i.e. the ratio between the peaks and the average value of the intensity in the ultrasound beam. A low BNR excludes the possibility of undesirably high energy concentrations in the beam.

ERA and BNR are measured according to the FDA methods. All treatment heads are provided with contact control.

TECHNICAL DATA

Mains voltage	:	115 V~ 50/60 Hz $\pm 15\%$ and 220-240 V~ 50/60 Hz $\pm 10\%$
Current consumption	:	max. 0.70 A at 115 V / 0.35 A at 230 V
Mains fuses	:	2 x T1.6H250V
Patient leakage current	:	typically 10 μ A (IEC requirement $\leq 100 \mu$ A)
Ditto, first fault condition	:	typically 30 μ A (IEC requirement $\leq 500 \mu$ A)
Earth (ground) leakage current	:	typically 60 μ A (IEC requirement $\leq 500 \mu$ A)
Ditto, single fault condition	:	typically 120 μ A (IEC requirement $\leq 1000 \mu$ A)
Earth resistance	:	typically 0.12 Ohm (IEC requirement ≤ 0.2 Ohm)
Medical device classification	:	IIb (according to Medical Device Directive (93/42/EEC))
Safety class	:	I* type BF**, according IEC 60601-1
Weight and dimensions	:	6.1 kg / 36 x 33.5 x 12 cm
Safety tests	:	TÜV Rheinland

Environment conditions for transport and storage

Environment temperature	:	-10° till +50° C
Relative humidity	:	10 till 95 %
Atmospheric pressure	:	500 till 1060 hPa

Environments conditions for normal use

Environment temperature	:	10° till 40° C
Relative humidity	:	10 till 90 %
Atmospheric pressure	:	500 till 1060 hPa

Medical class IIb

This equipment complies with all requirements of the Medical Device Directive (93/42/EEC)

International Safety class IEC60601-1

* I : The equipment has a safety earth (ground) connection, and must be connected to an earthed (grounded) wall socket.

** BF : The equipment has a floating patient circuit. The Sonopuls 590 meets the requirements of IEC 60601-1 and IEC 60601-2-5.

Copies of the test reports are available on request. Technical modifications reserved.

ORDERING DATA

For the ordering data of the Sonopuls, standard accessories and additional accessories we refer to the Catalogue for Physiotherapy.

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