



**Supply and Service of Medical Devices**

Unit 6B, Mulberry Trading Estate, Foundry Lane,  
Horsham, West Sussex, RH13 5PX.

Tel: 01403 597597 www.trimbio.co.uk

Email: serviceadmin@trimbio.co.uk

<b>Site:</b> Trimbio Workshop	<b>Device SKU:</b> SE10485
-------------------------------	----------------------------

**Secondhand Device Service Inspection Record**

<b>Schedule Date</b>	07/05/2026	<b>Scheduled Time</b>	07/05/2026 16:30
<b>Job No.</b>	J6322		

<b>Device</b>	<b>Serial No.</b>	<b>Manufacturer</b>	<b>Model</b>	<b>Model Number</b>
SE10485	78311470970	Laerdal	Suction Unit	

Electrical Safety Checks	Earth $\Omega$	Mains - PE M $\Omega$	Mains - AP M $\Omega$	Equipment Leakage		Applied Part Leakage	
				For $\mu$ A	Rev $\mu$ A	For $\mu$ A	Rev $\mu$ A
Medical-IEC62535	☐	N/A	N/A	22	26	5.4	6.4

<b>Medical- IEC62353</b>	<b>Portable Appliance Test</b>	<b>Function Test</b>	<b>Overall Pass or Fail</b>
Pass		Pass	Pass

<b>Device Notes</b>	
---------------------	--

<b>The above work has been carried out as per Trimbio's standard service terms and conditions and procedures by:</b>	<b>Engineers Name</b>	Ilias Savvinidis
--	-----------------------	------------------